

The „Triple A“ of Bone Regeneration: A for all indications, A for all users, A for all patients

Directly injectable from the syringe

Putty-like consistency

In situ hardening

Open interconnected porosity

100% synthetic

Osseo-regenerative

Large indication spectrum

Optimal bone regeneration



easy-graft™

Is a bioresorbable, completely synthetic bone graft for bone defects. It consists of two components: Granules and BioLinker™. When combining these two components, easy-graft™ becomes a putty-like biomaterial and can be applied directly from the syringe into the bone defect.



easy-graft™

In contact with body fluids – such as blood or saliva – the material hardens within minutes and forms a stable, porous bone substitute material. Due to the porosity of the material, the absorption of blood is possible and thus positively influences the healing process.



easy-graft™

Consists only of known and well-proven biomaterials: High-purity β -Tri-Calciumphosphate and Polylactide.



easy-graft™


The biomaterial resorbs completely in 9 to 15 months and will be replaced by autologous bone. easy-graft™ is a biocompatible and osseo-regenerativ material. The implantation will be proceed after 8 to 12 weeks.

Publications on our biomaterials: Thoma, K. Schmidlin, P. et al., Oral Surg Oral Med Oral Pathol Oral Radiol Endod, 2006, in press. Nair, PNR. et al., Journal of Biomaterials Applications, Volume 00 –2006, in press. Schmidlin, P. et al., Schweiz Monatsschr Zahnmed, Vol 114: 4 / 2004, p. 328 – 336. Schug, J. et al., Die Zahnarztwoche (DZW) 47 / 2002, p. 14-15. Nair, PNR. et al., Journal of Oral Surgery, Oral Medicine, Oral Pathology, 2004; Vol. 97: p. 559-569

Products	Granulate size	Indikation
easy-graft™ 150 6x 0,15 ml REF C11-013	500-600 μ m	Periodontal defects and small defects
easy-graft™ 400 3x 0,4ml REF C11-002	500-1000 μ m	Implatations, Alveolus, Augmentations, Cysts Sinus Lift
easy-graft™ 400 6x 0,4ml REF C11-003	500-1000 μ m	Implatations, Alveolus, Augmentations, Cysts Sinus Lifts



Degradable Solutions AG
Wagistrasse 23 · CH-8952 Schlieren ·
Tel: +41 (0)43 433 62 60 · Fax: +41 (0)43 433 62 61
E-Mail: kinzer@degradable.ch · www.degradable.com



DS
DENTAL

AAA

The „Triple A“
of Bone Regeneration



*easy-graft*TM

injectable

in situ hardening

100% synthetic bonegraft



sticky granules

bionic